

# Elastocon

## Accessories for the Elastocon Cell Oven



**Stand, EB 011**, stand to support the sample holder while mounting test pieces for ageing.

**Compression Set Rig, EV 03**, according to ISO 815. The Rig consists of two circular polished plates of stainless steel with a central screw and spacer. Designed to be suitable for holding in a conventional vice to aid sample loading, the rig utilises the minimum of components. Spacers for both small and large test pieces are included.

A spacer set with spacers from 1 to 20 mm in 0,5 mm step is also available.

**Tension Set Rig, EV 04**, according to ISO 2285.

The Rig is made of stainless steel and can be adjusted from 25 to 100% extension. The rig can easily be mounted on the sample holder of the Cell Oven.

**ARTEC**  
TESTNOLOGY

ARTEC TESTNOLOGY  
P.O. Box 12 - NL-5330 AA KERKDRIEL  
Salie 15 - NL-5331 DJ KERKDRIEL  
Tel: +31(0)418 637590 - Fax: +31(0)418 637599  
E-mail: [info@artec.nl](mailto:info@artec.nl) - Web: <http://www.artec.nl>